

Dr. PACKER
Phaco Chop
AMO Signature
February 2009

<i>Phaco tip: Gold 20 Ga 30 degree tip</i>	<i>Incision:</i>	<i>Disposable</i>
--	------------------	-------------------

↓ **Start here**

Memory	Phaco 1 CHOP	Phaco 2 SEGMENT	Phaco 3 EPINUCLEUS	IA Cortex Max	Viscoat Removal Co-Axial Curved Silicone Tip
Power	<i>50 Linear cont</i>	<i>100 Linear</i>	<i>20/20 Linear</i>		
Whitestar	<i>Long pulse 3 pps 6/12 (33%, 55 pps)</i>	<i>Long Pulse 6 pps</i>	<i>WS Variable Set 1</i>		
Ellips	<i>Off</i>	<i>off</i>	<i>On</i>		
Unocc/occ mode	<i>Off</i>	<i>off</i>	<i>On</i>		
Vacuum	<i>405 panel</i>	<i>350 linear</i>	<i>50)/300 Threshold/limit linear</i>	<i>500 linear</i>	<i>600 panel</i>
Flow	<i>36 panel</i>	<i>42 linear</i>	<i>34/26 linear</i>	<i>28 linear</i>	<i>50 linear</i>
Ramp	<i>40</i>	<i>40</i>	<i>40</i>	<i>60</i>	<i>60</i>
bottle ht (inches)	<i>106</i>	<i>106</i>	<i>105</i>	<i>106</i>	<i>106</i>

Vitrectomy

Hi-speed pneumatic	<i>Flow 20</i>	<i>Vacuum 250</i>	<i>Cut rate 1500</i>	<i>Bottle 40</i>
-------------------------------	----------------	-------------------	----------------------	------------------